

A Multiplex RT-LAMP-Lateral Flow Platform For Comprehensive Dengue Serotype Detection

Briana Gonzalez Bastida, Alina Darchieva, Kayla Jerome, Raniyah Anasatasia Le'iato, Aakriti Sapkota

Intent & Motivation

Objective: To develop a diagnostic tool that uses one-pot RT-LAMP setup for amplification and a multiplex LFA system to detect and differentiate all four dengue serotypes.

Motivation: To create a low-cost, reliable, and stable point-of-care diagnostic tool for dengue that allows early detection and serotype differentiation

Background

- Dengue fever is a **mosquito-borne viral disease** that affects **350 million people annually**, primarily in **tropical regions**, and is expanding globally due to rising temperatures and urbanization.¹
- Current diagnostics (RT-PCR, ELISA) require specialized lab equipment, **limiting accessibility in low-resource settings**.²
- **Reverse-transcription loop-mediated isothermal amplification (RT-LAMP)** enables rapid viral RNA amplification at a constant temperature.³
- **Lateral flow assays (LFAs)** provide a simple, paper-based visible line when the target molecule is detected.⁴
- Integrated RT-LAMP with gold nanoparticle (AuNP)-based LFAs offers a **sensitive, low-cost, point-of-care dengue diagnostics for under-resourced regions**.




Figure 1. Blood is collected and added to a tube with freeze-dried buffer, followed by RT-LAMP amplification at ~60–65 °C for 30–40 min. The amplified products are transferred to an LFA strip for serotype-specific visual detection, where positive results show test lines along with a control line, and negative results show only the control line.

References & Acknowledgements



QR code links to references

This research was supported by the VIP (Vertically Integrated Projects) Program at the University of Maryland. We would like to thank Dr. Spirito, peer mentors and the White Lab for providing resources and guidance throughout the project. Figures included in this work were created using Biorender.

Methods

Reverse-Transcription LAMP




Figure 2: RT-LAMP Process. Binding sites for the LAMP primers are highlighted in the figure.

- Target RNA is converted to cDNA and amplified by multiplex LAMP targeting all four dengue serotypes⁵
- LAMP utilizes 4-6 primers and Bst 2.0 polymerase
- Colorimetric and fluorescent LAMP were performed for initial testing

One Pot Heating System




Figure 3: 3D Printed One-Pot Chemical Heating system (1 square = 1 cm²)

- Custom 3D-printed acrylonitrile butadiene styrene (ABS) lid and base hold a test tube over a Mg-Fe alloy chamber that generates and maintains 60-65°C⁶ for LAMP upon water activation

Multiplex LFA Layout




Figure 4: Layout of the multiplex LFA for dengue serotype detection

- An amplified sample is applied to a pad that isolates LAMP products
- LAMP products bind to complementary capture probes immobilized on the test lines, that are designed by 2D modeling of LAMP products to ensure complementary binding
- Streptavidin-coated gold nanoparticles bind to biotin-tagged LAMP products to produce a visible signal
- A control probe is included in the control line to confirm proper function
- High salt concentrations enable probe binding without antibodies use⁷

Results & Discussion




Figure 5: Fluorescent LAMP ~ Optimization of LAMP amplification conditions for dengue detection using synthetic DNA.

- The new primer mix with betaine (**A**) demonstrates earlier amplification in the 10⁻² sample (0.0775 nM) compared to the old primer mix (**B**). The 10⁻⁴ sample (0.000775 nM) did not show clear amplification and was dominated by background noise. Negative controls exhibited minor increases in fluorescence, likely due to non-specific amplification or experimental variability.


Streptavidin AuNPs

Figure 6: Streptavidin Gold Nanoparticle (AuNP) Synthesis and Conjugation Optimization

- 12 nm AuNPs were synthesized via citrate reduction
- Conjugation was optimized by adjusting pH, BSA, and streptavidin concentration
- The reaction was scaled up 5x for downstream use



Pellet (before resuspension)



Conjugates (AuNP + Streptavidin)

Discussion:

- Successful fluorescence was seen in the 0.0775 nM sample with betaine added in under 40 minutes.
- Streptavidin AuNP synthesis was achieved at 10 mg/μL, with proteins remaining reddish before and after 10% NaCl addition, confirming successful conjugation

Future Work

- Evaluate one-pot heating system stability
- Validate LFA performance using salt-mediated immobilization
- Conduct cost analysis of system components
- Expand testing to all dengue serotypes beyond DENV-2